

Sample collection, immunostaining and optical clearing

 Hernán Morales-Navarrete  Yannis Kalaidzidis  Frank Jülicher  Benjamin M Friedrich  Marino Zerial

Updated date: Jan 2, 2020

 An abbreviated version of this protocol was published in eLIFE in Jun 2019

Liquid-crystal organization of liver tissue

DOI: 10.7554/eLife.44860

Detailed protocol

Fluorescence immunostaining - floating sections

A) Vibratome

Mold

1. Turn on the water bath at 60°C.
2. Prepare 4% agarose in PBS and maintain it melted in the bath water at 60°C. Put agarose in the plastic mold and add one piece of the dry liver at the bottom.
3. Make vibratome sections (100 µm thickness). Add PBS 500 µl/well in a 48 well plate. Just 1 slice per well.

B) Immunofluorescence

1. Remove the agarose from the tissue and permeabilize with 0.5% Triton X-100 in PBS (60 minutes) (300 µl/well).
2. Add primary antibody in TxBuffer (2 overnights room temperature)
3. Wash with 0.3% Triton/PBS (5 times 15 minutes)
4. Add secondary antibody + DAPI in TxBuffer (2 overnights room temperature)
5. Wash with 0.3% Triton/PBS (5 times 15 minutes)
6. Wash in PBS (3 times 1 minute)

C) M-SeeD Clearing.

1. Add 200 µl of 25% fructose for 4 hrs,
2. then 50% fructose for 4 hrs ,
3. 75% fructose ON
4. 100% fructose ON.
5. Add 200 µl of SeeD ON. All these steps are at room temperature.
6. Mount on a glass slide with SeeD solution. #1.5 coverslips (thickness 0.17 ± 0.005 mm). RI:1,49

- Different concentrations of fructose are prepared diluting 100% fructose with water.

Buffer recipes

TxBuffer (1L):

0.2% gelatin

300mM NaCl

0.3% Triton X-100

Add PBS to make 1L, aliquot to 50mL falcon tubes, store at -20°C

Heat up PBS, gelatin and NaCl to dissolve. After that when the solution is under 40°C add Triton.

Modified SeeDB clearing

Reagents

- M-SeeDB

80.2% (wt/wt) fructose, 0.5% 1-thioglycerol, ~0.1M phosphate buffer (pH7.5)*

- 100% Fructose pH7.5

100% (wt/v) fructose, 0.5% 1-thioglycerol, 0.1M phosphate buffer (pH7.5)

M-SeeDB recipe

Fructose 40.1 g

1M Na₂HPO₄ 2.4 ml

1M NaH₂PO₄ 0.6 ml

1-thioglycerol 0.15 ml

ddH₂O up to 50 g

100% Fructose

Fructose 40 g

1M Na₂HPO₄ 3.2 ml

1M NaH₂PO₄ 0.8 ml

1-thioglycerol 0.2 ml

ddH₂O up to 40 mL

* To keep fluorescent signal, the solution has to be buffered.

How to cite: (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Morales-Navarrete, H. , Kalaidzidis, Y. , Jülicher, F. , Friedrich, B. M. and Zerial, M. (2020). Sample collection, immunostaining and optical clearing.

Bio-protocol Preprint. [bio-protocol.org/prep173](https://doi.org/10.21203/rs.3.rs-173).

2. Morales-Navarrete, H., Nonaka, H., Scholich, A., Segovia-Miranda, F., de Back, W., Meyer, K., Bogorad, R. L., Koteliensky, V., Brusch, L., Kalaidzidis, Y., Jülicher, F., Friedrich, B. M. and Zerial, M.(2019). Liquid-crystal organization of liver tissue. eLIFE. DOI: [10.7554/eLife.44860](https://doi.org/10.7554/eLife.44860)

Copyright: Content may be subjected to copyright.